


E 989

Universal bracket

for 8 insoles for theatre shoes, for example.



- For use in insert E 984

EAN: 4002512431828 / Material number: 05376740 / Old Material Number: 69598901D

Product category	
Large-capacity washer-disinfectors, medical	•
Machine affiliation	
PG 8528	•
Product category	
Other inserts	•
Product characteristics	
Material	Stainless steel
Colour	Stainless Steel
Application area	
Reprocessing of OR shoes	•
Capacity	
Insoles [number]	8
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in mm	0
External dimensions, net width in mm	0
External dimensions, net depth in mm	0
External dimensions, gross height in mm	0 7/8 (20)
External dimensions, gross width in mm	13 (330)
External dimensions, gross depth in in (mm)	16 3/8 (415)
Net weight in kg	1 (0.10)
Gross weight in kg	1 (0.26)