


E 747

Module

for optimum loading of centrifuge tubes, vials and test tubes, for example.

EAN: 4002512591218 / Material number: 05464630 / Old Material Number: 69574701D



- 104 nozzles for centrifuge tubes, phials, test tubes, etc.
- 104 nozzles 2.5 x 110 mm

Product category	
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
G 7825	•
Product category	
Insert / module for test tubes	•
Product characteristics	
Material	Stainless steel, Plastic
Colour	Stainless Steel, Black
Application area	
Reprocessing of pipettes / phials	•
Capacity	
Test tubes [number]	104
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in in. (mm)	6 5/8 (168)
External dimensions, net width in in. (mm)	19 3/8 (492)
External dimensions, net depth in in. (mm)	19 9/16 (496)
External dimensions, gross height in in. (mm)	25 5/8 (650)
External dimensions, gross width in in. (mm)	26 (660)
External dimensions, gross depth in in. (mm)	31 15/16 (810)
Net weight in lb (kg)	84 (38.00)
Gross weight in lb (kg)	98 (44.00)