


E 746
Module
for the optimum loading of up to 23 pipettes.



- Suitable for 23 pipettes arranged diagonally
- 10 pipettes up to 560 mm and 13 pipettes up to 490 mm
- Distance between holders 20 or 26 mm

EAN: 4002512591195 / Material number: 05555280 / Old Material Number: 69574601D

Product category	
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
G 7825	•
Product category	
Insert / module for pipettes	•
Product characteristics	
Material	Stainless steel, Plastic
Colour	Stainless Steel, Black
Capacity	
Pipettes [number]	23
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in in. (mm)	13 (330)
External dimensions, net width in in. (mm)	19 3/8 (492)
External dimensions, net depth in in. (mm)	19 9/16 (496)
External dimensions, gross height in in. (mm)	15 3/4 (400)
External dimensions, gross width in in. (mm)	23 1/4 (590)
External dimensions, gross depth in in. (mm)	21 5/16 (540)
Net weight in lb (kg)	55 (24.60)
Gross weight in lb (kg)	56 (25.00)