

E 469

Irrigation sleeves

for the optimum loading of MIS instruments from urology, L 300 mm, Ø 11 mm.



- Cap, opening Ø 6 mm
- With mounting clip

EAN: 4002512598446 / Material number: 05563950 / Old Material Number: 69546901D

Product category	
Washer-disinfectors, medical	•
Large-capacity washer-disinfectors, medical	•
Machine affiliation	
G 7823	•
G 7824	•
G 7825	•
G 7826	•
PG 8527	•
PG 8528	•
PG 8535	•
PG 8536	•
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8592	•
Product category	
Irrigation sleeve	•
Product characteristics	
Material	Stainless steel, Plastic
Colour	Stainless Steel, Black
Application area	
Reprocessing of MIS instruments	•
Reprocessing of lumened instruments	•
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in mm	300
External dimensions, net width in mm	0
External dimensions, net depth in mm	0
External dimensions, gross height in mm	5 3/4 (145)
External dimensions, gross width in mm	10 1/2 (265)
External dimensions, gross depth in in (mm)	6 1/8 (155)
Diameter in in (mm)	0 1/4 (6)
Net weight in kg	1 (0.05)
Gross weight in kg	1 (0.05)