


E 523

Insert

for the optimum loading of 6 mesh trays.



- Equipped with 7 holders (6 compartments), H 145, spacing 50 mm

EAN: 4002513199116 / Material number: 05649700 / Old Material Number: 69552301D

<b>Product category</b>	
Thermal disinfectors	•
<b>Machine affiliation</b>	
PG 8581	•
PG 8591	•
<b>Product category</b>	
Insert for trays	•
<b>Product characteristics</b>	
Material	Stainless steel
Colour	Stainless Steel
Number of compartments	6
<b>Application area</b>	
Reprocessing of dental instruments	•
<b>Standard electrical connection</b>	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
<b>Dimensions and weight</b>	
External dimensions, net height in mm	151
External dimensions, net width in mm	220
External dimensions, net depth in mm	450
External dimensions, gross height in mm	7 3/4 (195)
External dimensions, gross width in mm	9 1/2 (240)
External dimensions, gross depth in in (mm)	19 1/2 (495)
Net weight in kg	4 (1.40)
Gross weight in kg	4 (1.65)