

E 801/1

Insert

for the optimum loading of 8 mouthwash glasses.



- Suitable for 8 mouth rinse tumblers, max. Ø 75 mm
- Plastic-coated

EAN: 4002513220650 / Material number: 05673820 / Old Material Number: 68580102D

Product category	
Thermal disinfectors	•
Machine affiliation	
G 7831	•
PG 8535	•
PG 8536	•
PG 8581	•
PG 8591	•
PWD 8531	•
PWD 8532	•
Product category	
Other inserts	•
Product characteristics	
Material	Wire, rilsan-coated
Colour	White
Application area	
Reprocessing of dental instruments	•
Capacity	
Mouth rinse cups [number]	8
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in mm	155
External dimensions, net width in mm	200
External dimensions, net depth in mm	320
External dimensions, gross height in mm	10 3/4 (270)
External dimensions, gross width in mm	20 1/8 (510)
External dimensions, gross depth in in (mm)	9 1/2 (240)
Net weight in kg	3 (1.00)
Gross weight in kg	4 (1.50)