

E 476

Holder

for micro instruments with a Ø of 4-8 mm.



- For use in trays with 5 mm mesh
- 50 items per bag

EAN: 4002513239454 / Material number: 05692400 / Old Material Number: 69747601D

Product category	
Thermal disinfectors	•
Washer-disinfectors, medical	•
Large-capacity washer-disinfectors, medical	•
Machine affiliation	
G 7823	•
G 7824	•
G 7825	•
G 7826	•
PG 8527	•
PG 8528	•
G 7831	•
PG 8535	•
PG 8536	•
PG 8581	•
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8591	•
PG 8592	•
PWD 8531	•
PWD 8532	•
Product category	
Holder	•
Product characteristics	
Material	Plastic
Colour	Grey
Application area	
Reprocessing of MIS instruments	•
Reprocessing of ophthalmology instruments	•
Reprocessing of dental instruments	•
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in mm	0
External dimensions, net width in mm	0
External dimensions, net depth in mm	0
External dimensions, gross height in mm	0 1/2 (12)
External dimensions, gross width in mm	3 3/16 (80)
External dimensions, gross depth in in (mm)	4 3/4 (120)
Net weight in kg	1 (0.07)
Gross weight in kg	1 (0.08)