


A 19

Lid

For fixing delicate wash items.



- Perforated plate 7 x 7 x 3 mm
- Suitable for insert E 969

EAN: 4002513298758 / Material number: 05746210 / Old Material Number: 69796901D

Product category	
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PG 8527	•
PLW 8615	•
PLW 8616	•
PLW 8617	•
Product category	
Lid	•
Product characteristics	
Material	Stainless steel
Colour	Stainless Steel
Application area	
Reprocessing of laboratory glassware	•
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in mm	18
External dimensions, net width in mm	251
External dimensions, net depth in mm	351
External dimensions, gross height in mm	2 7/8 (70)
External dimensions, gross width in mm	10 7/8 (275)
External dimensions, gross depth in in (mm)	20 3/8 (515)
Net weight in kg	2 (0.55)
Gross weight in kg	2 (0.70)