

E 131/1

Insert

for optimum loading of 5 mesh and kidney trays.



- Suitable for 5 mesh or kidney trays
- Equipped with 6 holders, height: 6 3/8 (160) in (mm) mm, distance: 3 1/4 (80) in (mm) mm

EAN: 4002513317336 / Material number: 05771030 / Old Material Number: 69513102D

<b>Product category</b>	
Washer-disinfectors, medical	•
<b>Machine affiliation</b>	
PG 8535	•
PG 8536	•
PG 8581	•
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8591	•
PG 8592	•
<b>Product category</b>	
Insert for mesh trays	•
Insert for kidney dishes / basins	•
<b>Product characteristics</b>	
Material	Stainless steel
Colour	Stainless Steel
Distance between rails in mm	3 1/4 (80)
Rail height in mm	6 3/8 (160)
<b>Application area</b>	
Reprocessing of kidney dishes	•
Reprocessing of surgical instruments	•
<b>Capacity</b>	
Kidney dishes [number]	5
<b>Standard electrical connection</b>	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
<b>Dimensions and weight</b>	
External dimensions, net height in mm	168
External dimensions, net width in mm	180
External dimensions, net depth in mm	495
External dimensions, gross height in mm	7 3/8 (185)
External dimensions, gross width in mm	7 7/8 (200)
External dimensions, gross depth in in (mm)	20 7/8 (530)
Net weight in kg	2 (0.61)
Gross weight in kg	2 (0.90)