

E 492

Insert

for the optimum loading of kidney trays.



- equipped with 9 holders, height 3 1/2 (86) in (mm) mm, distance 2 (49) in (mm) mm
- For use in lower basket

EAN: 4002513317343 / Material number: 05771140 / Old Material Number: 69549201D

Product category	
Washer-disinfectors, medical	•
Machine affiliation	
G 7831	•
PG 8535	•
PG 8536	•
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8592	•
Product category	
Insert for kidney dishes / basins	•
Other inserts	•
Product characteristics	
Material	Stainless steel
Colour	Stainless Steel
Distance between rails in mm	2 (49)
Rail height in mm	3 1/2 (86)
Application area	
Reprocessing of kidney dishes	•
Capacity	
Kidney dishes [number]	9
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in mm	120
External dimensions, net width in mm	256
External dimensions, net depth in mm	474
External dimensions, gross height in mm	5 (125)
External dimensions, gross width in mm	12 5/8 (320)
External dimensions, gross depth in in (mm)	19 (480)
Net weight in kg	3 (1.06)
Gross weight in kg	3 (1.31)