


E 965

Insert

for the optimum loading of beakers 250 - 600 ml.



- Equipped with 15 x 3 stay bars for beakers from 250 ml to 600 ml

EAN: 4002513369229 / Material number: 05848290 / Old Material Number: 69596501D

Product category	
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PG 8527	•
Product category	
Insert / module for beakers	•
Product characteristics	
Material	Stainless steel
Colour	Stainless Steel
Application area	
Reprocessing of laboratory glassware	•
Capacity	
Beakers [quantity]	98
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in mm	173
External dimensions, net width in mm	357
External dimensions, net depth in mm	522
External dimensions, gross height in mm	11 1/16 (280)
External dimensions, gross width in mm	21 5/16 (540)
External dimensions, gross depth in in (mm)	22 1/16 (560)
Net weight in kg	5 (2.07)
Gross weight in kg	7 (2.86)