

U 168/1

Lower basket

for the ideal loading of theatre shoes.



- Suitable for 20 theatre shoes up to size 45
- 20 holders 11 5/8 (295) in (mm) mm
- Bracket for ML/2 for wash cart recognition / spray arm control

EAN: 4002513390940 / Material number: 05878590 / Old Material Number: 69116803D

Product category	
Washer-disinfectors, medical	•
Machine affiliation	
PG 8535	•
PG 8536	•
Product category	
Insert for OR shoes	•
Lower basket	•
Product characteristics	
Material	Stainless steel
Colour	Stainless Steel
Rail height in mm	11 5/8 (295)
Application area	
Reprocessing of OR shoes	•
Capacity	
OR shoes [number]	20
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in in. (mm)	12 7/16 (315)
External dimensions, net width in in. (mm)	21 1/8 (535)
External dimensions, net depth in in. (mm)	20 5/16 (515)
External dimensions, gross height in in. (mm)	12 5/8 (320)
External dimensions, gross width in in. (mm)	21 1/2 (545)
External dimensions, gross depth in in (mm)	21 1/2 (545)
Net weight in lb (kg)	8 (3.20)
Gross weight in lb (kg)	12 (5.20)