

O 177/1

Upper basket

for optimum cleaning with 26 nozzles and holding of various inserts.



- Right side for different inserts
- Left side with 26 nozzles and 7 funnels supplied loose
- Equipped with height adjustable frame
- Vertical clearance 230/205 mm
- Height adjustable: – 20 / – 40 mm

EAN: 4002513412413 / Material number: 05913050 / Old Material Number: 69117702D

Product category	
Thermal disinfectors	•
Washer-disinfectors, medical	•
Machine affiliation	
PG 8535	•
PG 8536	•
Product category	
Carrier / basket carrier	•
Upper basket	•
Product characteristics	
Material	Stainless steel, Plastic
Colour	Stainless Steel, Black
Application area	
Reprocessing of ENT instruments	•
Reprocessing of dental instruments	•
Capacity	
Handpieces [number]	11
Lumened instruments [number]	22
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in in. (mm)	10 3/8 (263)
External dimensions, net width in in. (mm)	19 5/8 (498)
External dimensions, net depth in in. (mm)	17 15/16 (455)
External dimensions, gross height in in. (mm)	11 1/16 (280)
External dimensions, gross width in in. (mm)	21 1/2 (545)
External dimensions, gross depth in in (mm)	22 1/16 (560)
Net weight in lb (kg)	9 (3.70)
Gross weight in lb (kg)	10 (4.52)