

E 441/1

Insert

for the optimum loading of Micro instruments.



- The mesh width on the base is 1.7 mm
- The sides are closed, stackable
- Sections w/ 6 adjustable dividers for positioning of instruments
- Can be steam sterilized at 121 °C/134 °C

EAN: 4002513460476 / Material number: 05988240 / Old Material Number: 69544102D

Product category	
Thermal disinfectors	•
Washer-disinfectors, medical	•
Machine affiliation	
G 7831	•
PG 8535	•
PG 8536	•
PG 8581	•
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8591	•
PG 8592	•
Product category	
Insert / module for MIS instruments	•
Product characteristics	
Material	Stainless steel, Plastic
Colour	Stainless Steel, Black
Application area	
Reprocessing of MIS instruments	•
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in mm	2 3/8 (60) in (mm)
External dimensions, net width in mm	7 1/4 (183) in (mm)
External dimensions, net depth in mm	11 3/16 (284) in (mm)
External dimensions, gross height in mm	3 15/16 (100)
External dimensions, gross width in mm	8 5/16 (210)
External dimensions, gross depth in in (mm)	12 7/16 (315)
Net weight in kg	2 (0.77)
Gross weight in kg	3 (0.92)