

FP 20

Filter plates

for the reprocessing of contra angles and hand pieces when using AUF 1/AUF 2.

EAN: 4002513505498 / Material number: 06059480 / Old Material Number: 69751201D



- To avoid blockages by micro particles
- Porosity 2
- 20 filter plates per pack

Product category	
Thermal disinfectors	•
Large-capacity washer-disinfectors, medical	•
Machine affiliation	
G 7823	•
G 7824	•
G 7825	•
G 7826	•
PG 8527	•
PG 8528	•
G 7831	•
PG 8535	•
PG 8536	•
Product characteristics	
Colour	White
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in mm	1 3/16 (30) in (mm)
External dimensions, net width in mm	3 15/16 (100) in (mm)
External dimensions, net depth in mm	4 15/16 (125) in (mm)
External dimensions, gross height in mm	1 1/8 (28)
External dimensions, gross width in mm	4 15/16 (125)
External dimensions, gross depth in in (mm)	6 5/8 (167)
Diameter in in (mm)	1 3/16 (30)
Net weight in kg	1 (0.06)
Gross weight in kg	1 (0.10)