


E 902/1

Mobile unit

for holding up to 3 MIS theatre sets in a modular system.



- Suitable for 6 insert modules with integrated nozzles/adapters
- Upper & lower levels free for inserts, e.g. E 142
- Equipped with connection for hot air drying
- Magnetic strip for wash cart recognition (incl. magnets ML)
- Simultaneous cleaning of mesh trays and MIS theatre sets

EAN: 4002513505542 / Material number: 06059530 / Old Material Number: 69590202D

Product category	
Large-capacity washer-disinfectors, medical	•
Machine affiliation	
PG 8527	•
PG 8528	•
Product category	
Mobile unit for MIS instruments	•
Product characteristics	
Material	Stainless steel
Colour	Stainless Steel
Application area	
Reprocessing of MIS instruments	•
Capacity	
DIN mesh trays [number]	3
SPRI mesh trays [number]	1
MIS sets [number]	3
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in in. (mm)	24 7/16 (620)
External dimensions, net width in in. (mm)	25 1/4 (640)
External dimensions, net depth in in. (mm)	31 1/8 (790)
External dimensions, gross height in in. (mm)	25 5/8 (650)
External dimensions, gross width in in. (mm)	25 5/8 (650)
External dimensions, gross depth in in (mm)	31 1/2 (800)
Net weight in lb (kg)	53 (24.00)
Gross weight in lb (kg)	75 (34.00)
Accessory pack	
E 362 blind stopper [number]	6
E 444 drum for fibre optic cables / suction hoses [number]	3
E 445 – 12 irrigation sleeve covers E 442, opening with Ø 6 mm [number]	2
E 446 – 12 irrigation sleeve covers E 443, opening with Ø 10 mm [number]	1
E 905/1 module insert for short MIS instruments [number]	3
E 906/1 module insert for long MIS instruments [number]	3
2 m silicone hose, Ø 5 mm	•
Magnetic strip for automatic program recognition (incl. magnets) [number]	1