


E 745/1

Module

for the optimum loading of up to 104 pipettes.



- Suitable for 104 pipettes up to 540 mm
- Holding frame compartment size 16 x 16 mm

EAN: 4002513735406 / Material number: 06233580 / Old Material Number: 69574502D

Product category	
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
G 7825	•
Product category	
Insert / module for pipettes	•
Product characteristics	
Material	Stainless steel, Plastic
Colour	Stainless Steel, Black
Application area	
Reprocessing of pipettes / phials	•
Capacity	
Pipettes [number]	104
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in in. (mm)	11 3/8 (288)
External dimensions, net width in in. (mm)	19 3/8 (492)
External dimensions, net depth in in. (mm)	19 9/16 (496)
External dimensions, gross height in in. (mm)	16 3/8 (415)
External dimensions, gross width in in. (mm)	20 7/8 (530)
External dimensions, gross depth in in (mm)	23 1/4 (590)
Net weight in lb (kg)	55 (24.60)
Gross weight in lb (kg)	56 (25.00)