

E 529

Mobile unit

for optimum loading of 4-5 ophthalmic theatre sets.



- 2 levels for mesh trays
- 1 level with 16 or 30 injector nozzles
- Connection option for E 478/1 to the inlet pipe
- Equipped with connection for hot air drying
- Magnetic strip for mobile unit recognition (without magnets ML)

EAN: 4002513792553 / Material number: 06415950 / Old Material Number: 69552901D

| | |
|---|-----------------|
| Product category | |
| Large-capacity washer-disinfectors, medical | • |
| Machine affiliation | |
| G 7823 | • |
| G 7824 | • |
| Product category | |
| Mobile unit for ophthalmology instruments | • |
| Product characteristics | |
| Material | Stainless steel |
| Colour | Stainless Steel |
| Application area | |
| Reprocessing of ophthalmology instruments | • |
| Capacity | |
| Surgery sets for micro instruments [number] | 4-5 |
| Lumened instruments [number] | 46 |
| Standard electrical connection | |
| Number of phases | 0 |
| Voltage in V | 0 |
| Frequency in Hz | 0 |
| Total rated load in kW | 0.00 |
| Fuse rating in A | 0 |
| Dimensions and weight | |
| External dimensions, net height in in. (mm) | 18 9/16 (471) |
| External dimensions, net width in in. (mm) | 20 7/8 (530) |
| External dimensions, net depth in in. (mm) | 25 7/16 (645) |
| External dimensions, gross height in in. (mm) | 21 11/16 (550) |
| External dimensions, gross width in in. (mm) | 25 1/4 (640) |
| External dimensions, gross depth in in (mm) | 25 5/8 (650) |
| Net weight in lb (kg) | 31 (14.00) |
| Gross weight in lb (kg) | 47 (21.00) |