


U 175/1

Lower basket

for the optimum loading of narrow necked glassware.



- Equipped with 33 injector nozzles
- Vertical clearance: 6 3/4 (170) in (mm) mm*
- Only suitable for use with upper baskets O 175 / O 187 / O 184

* Depending on the laboratory glassware the max. vertical clearance can be lower.

EAN: 4002513805543 / Material number: 06439170 / Old Material Number: 69117503D

Product category	
Glassware Washers, laboratory	•
Machine affiliation	
PG 8536	•
Product category	
Injector basket	•
Lower basket	•
Product characteristics	
Material	Stainless steel
Colour	Stainless Steel
Maximum vertical clearance in mm	6 3/4 (170)
Application area	
Reprocessing of laboratory glassware	•
Capacity	
Narrow necked glassware [number]	33
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in in. (mm)	9 7/8 (250)
External dimensions, net width in in. (mm)	18 3/4 (475)
External dimensions, net depth in in. (mm)	20 15/16 (531)
External dimensions, gross height in in. (mm)	10 11/16 (270)
External dimensions, gross width in in. (mm)	20 11/16 (525)
External dimensions, gross depth in in (mm)	21 5/16 (540)
Net weight in lb (kg)	15 (6.40)
Gross weight in lb (kg)	17 (7.65)
Accessory pack	
E 351 injector nozzle, Ø 4 mm, length 160 mm [number]	33
E 353 clip for injector nozzle E 351 [number]	33
Rinse nozzle for powder dispenser [number]	1