


U 184/1

Lower basket

for the optimum loading of centrifuge tubes, vials and test tubes.



- Equipped with 96 injector nozzles
- Only suitable for use with upper baskets O 175 / O 187 / O 184
- Vertical clearance: Approx. 6 3/4 (170) in (mm) mm

EAN: 4002513805659 / Material number: 06439180 / Old Material Number: 69118403D

<b>Product category</b>	
Glassware Washers, laboratory	•
<b>Machine affiliation</b>	
PG 8536	•
<b>Product category</b>	
Injector basket	•
Lower basket	•
<b>Product characteristics</b>	
Material	Stainless steel, Plastic
Colour	Stainless Steel, Black
Maximum vertical clearance in mm	6 3/4 (170)
<b>Application area</b>	
Reprocessing of pipettes / phials	•
Reprocessing of test tubes	•
<b>Standard electrical connection</b>	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
<b>Dimensions and weight</b>	
External dimensions, net height in in. (mm)	9 7/8 (250)
External dimensions, net width in in. (mm)	20 5/16 (515)
External dimensions, net depth in in. (mm)	20 15/16 (531)
External dimensions, gross height in in. (mm)	10 11/16 (270)
External dimensions, gross width in in. (mm)	20 11/16 (525)
External dimensions, gross depth in in (mm)	21 5/16 (540)
Net weight in lb (kg)	12 (5.00)
Gross weight in lb (kg)	14 (6.00)
<b>Accessory pack</b>	
ID 90, injector nozzle, Ø 2.5 mm, length 90 mm [number]	96