

E 436/3

Mobile unit

for optimum loading of approx. 3 anaesthesia sets.



- Suitable for approx. 3 anaesthesia sets.
- 8 nozzles with spring holders for breathing tubes up to 1.5 m
- Storage on spiral
- Equipped with connection for hot air drying
- Bracket for ML/2 for wash cart recognition / spray arm control

EAN: 4002513805673 / Material number: 06439240 / Old Material Number: 69543604D

Product category	
Washer-disinfectors, medical	•
Machine affiliation	
PG 8536	•
Product category	
Mobile unit for anaesthesia instruments	•
Product characteristics	
Material	Stainless steel
Colour	Stainless Steel
Application area	
Reprocessing anaesthetic instruments	•
Capacity	
AN sets [number]	3
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in in. (mm)	20 (507)
External dimensions, net width in in. (mm)	21 1/8 (535)
External dimensions, net depth in in. (mm)	21 1/8 (535)
External dimensions, gross height in in. (mm)	20 1/8 (510)
External dimensions, gross width in in. (mm)	21 5/16 (540)
External dimensions, gross depth in in (mm)	20 1/2 (520)
Net weight in lb (kg)	15 (6.60)
Gross weight in lb (kg)	16 (7.00)
Accessory pack	
E 430/1 1/3 mesh tray [number]	1
E 431 injector nozzle for concertina hose, Ø 8 x 193 mm [number]	1
E 432 holder for concertina breathing tubes [number]	1
E 433 holder for silicone breathing tubes [number]	2
E 434 holder for paediatric breathing tubes [number]	1
E 466 injector nozzle 8 x 333 mm, for breathing bags [number]	6
E 496 injector nozzle 4 x 120 mm [number]	10