

E 427

Insert

for the optimum loading of 6 laryngoscopes.



- Suitable for 6 laryngoscopes

EAN: 4002513808209 / Material number: 06442120 / Old Material Number: 69542701D

<b>Product category</b>	
Washer-disinfectors, medical	•
Large-capacity washer-disinfectors, medical	•
<b>Machine affiliation</b>	
G 7823	•
G 7824	•
G 7825	•
G 7826	•
PG 8528	•
PG 8535	•
PG 8536	•
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8592	•
<b>Product category</b>	
Insert / module for anaesthetic instruments	•
<b>Product characteristics</b>	
Material	Stainless steel
Colour	Stainless Steel
<b>Application area</b>	
Reprocessing anaesthetic instruments	•
<b>Standard electrical connection</b>	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
<b>Dimensions and weight</b>	
External dimensions, net height in mm	92
External dimensions, net width in mm	210
External dimensions, net depth in mm	134
External dimensions, gross height in mm	6 1/2 (164)
External dimensions, gross width in mm	9 3/4 (245)
External dimensions, gross depth in in (mm)	9 7/8 (250)
Net weight in kg	2 (0.90)
Gross weight in kg	3 (1.00)