


U 874/1

Lower basket

for the optimum loading of different inserts and mesh trays.



- Vertical clearance can vary with different upper baskets
- With bracket for magnetic strip for ML/2 autom. cart recognition

EAN: 4002513808919 / Material number: 06442940 / Old Material Number: 69187403D

<b>Product category</b>	
Glassware Washers, laboratory	•
Washer-disinfectors, medical	•
<b>Machine affiliation</b>	
PG 8535	•
PG 8536	•
<b>Product category</b>	
Carrier / basket carrier	•
Lower basket	•
<b>Product characteristics</b>	
Material	Stainless steel
Colour	Stainless Steel
<b>Application area</b>	
Reprocessing of laboratory glassware	•
Reprocessing of MIS instruments	•
Reprocessing of GYN instruments	•
Reprocessing of ENT instruments	•
Reprocessing of OR shoes	•
Reprocessing of basins	•
Reprocessing of kidney dishes	•
Reprocessing of surgical instruments	•
Reprocessing of mesh trays	•
Reprocessing of different utensils	•
<b>Standard electrical connection</b>	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
<b>Dimensions and weight</b>	
External dimensions, net height in in. (mm)	2 (50)
External dimensions, net width in in. (mm)	20 15/16 (531)
External dimensions, net depth in in. (mm)	20 5/16 (515)
External dimensions, gross height in in. (mm)	3 15/16 (100)
External dimensions, gross width in in. (mm)	20 7/8 (530)
External dimensions, gross depth in in. (mm)	21 1/2 (545)
Net weight in lb (kg)	5 (2.10)
Gross weight in lb (kg)	7 (2.97)