

E 715/1

Mobile unit

for holding up to 4 AN sets and/or 12 breathing tubes (1.5 m) and accessories.



- Holds breathing tubes on a spiral rack
- Equipped with connection for hot air drying
- Magnetic strip for mobile unit recognition (without magnets ML)

EAN: 4002513872699 / Material number: 06491890 / Old Material Number: 69571502D

<b>Product category</b>	
Large-capacity washer-disinfectors, medical	•
<b>Machine affiliation</b>	
G 7825	•
G 7826	•
<b>Product category</b>	
Mobile unit for anaesthesia instruments	•
<b>Product characteristics</b>	
Material	Stainless steel, Plastic
Colour	Stainless Steel, Black
<b>Application area</b>	
Reprocessing anaesthetic instruments	•
<b>Capacity</b>	
AN sets [number]	4
<b>Standard electrical connection</b>	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
<b>Dimensions and weight</b>	
External dimensions, net height in in. (mm)	20 7/8 (530)
External dimensions, net width in in. (mm)	23 5/8 (600)
External dimensions, net depth in in. (mm)	25 7/16 (645)
External dimensions, gross height in in. (mm)	22 1/16 (560)
External dimensions, gross width in in. (mm)	25 1/4 (640)
External dimensions, gross depth in in (mm)	26 13/16 (680)
Net weight in lb (kg)	22 (9.70)
Gross weight in lb (kg)	23 (10.00)
<b>Accessory pack</b>	
E 430/1 1/3 mesh tray [number]	1
E 431 injector nozzle for concertina hose, Ø 8 x 193 mm [number]	1
E 432 holder for concertina breathing tubes [number]	1
E 433 holder for silicone breathing tubes [number]	3
E 434 holder for paediatric breathing tubes [number]	1
E 466 injector nozzle 8 x 333 mm, for breathing bags [number]	6
E 496 injector nozzle 4 x 120 mm [number]	15
Magnetic strip for automatic program recognition (without magnets) [number]	1