

E 479

Holder

for micro instruments up to 4 mm Ø.



- For use in trays with 5 mm mesh
- 50 items per bag

EAN: 4002513918175 / Material number: 06553750 / Old Material Number: 69747901D

Product category	
Washer-disinfectors, medical	•
Large-capacity washer-disinfectors, medical	•
Machine affiliation	
G 7823	•
G 7824	•
G 7825	•
G 7826	•
PG 8527	•
PG 8528	•
PG 8535	•
PG 8536	•
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8592	•
Product category	
Holder	•
Product characteristics	
Material	Plastic
Colour	Transparent
Application area	
Reprocessing of MIS instruments	•
Reprocessing of ophthalmology instruments	•
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in in. (mm)	0 1/4 (5)
External dimensions, net width in in. (mm)	0 7/8 (21)
External dimensions, net depth in in. (mm)	0 7/8 (21)
External dimensions, gross height in in. (mm)	3 3/4 (95)
External dimensions, gross width in in. (mm)	4 3/8 (110)
External dimensions, gross depth in in. (mm)	3 3/16 (80)
Net weight in lb (kg)	1 (0.19)
Gross weight in lb (kg)	1 (0.23)