

E 494
 Insert
 for the gentle holding of 5 microtiter plates.



EAN: 4002513926064 / Material number: 06570920 / Old Material Number: 69549401D

Product category	
Glassware Washers, laboratory	•
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
G 7825	•
PG 8527	•
PG 8535	•
PG 8536	•
PG 8504	•
PG 8583	•
PG 8583 CD	•
PG 8593	•
PLW 8615	•
PLW 8616	•
PLW 8617	•
PLW 6011	•
PLW 6111	•
PLW 8505	•
Product category	
Other inserts	•
Product characteristics	
Material	Stainless steel
Colour	Stainless Steel
Application area	
Reprocessing of laboratory glassware	•
Reprocessing of different utensils	•
Capacity	
Microtiter plates [number]	5
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in in. (mm)	1 7/16 (35)
External dimensions, net width in in. (mm)	8 1/8 (205)
External dimensions, net depth in in. (mm)	17 3/8 (440)
External dimensions, gross height in in. (mm)	3 15/16 (100)
External dimensions, gross width in in. (mm)	10 7/16 (265)
External dimensions, gross depth in in (mm)	7 1/8 (180)
Net weight in lb (kg)	2 (0.85)
Gross weight in lb (kg)	3 (1.00)