

E 405/1

Mobile unit

for optimum loading of up to 38 pipettes, connection for hot air drying.

EAN: 4002514087078 / Material number: 06727920 / Old Material Number: 69540502D



- Suitable for 38 pipettes in 3 rows
- 1. Row: 10 pipettes x 100 ml (up to 550 mm), bar distance 20 mm
- 2nd Row: 14 pipettes x 25 ml, bar distance 26 mm
- 3. Row: 14 pipettes x 10 ml, bar distance 26 mm
- Bracket for ML/2 for wash cart recognition / spray arm control

Product category	
Glassware Washers, laboratory	•
Machine affiliation	
PG 8535	•
PG 8536	•
Product category	
Mobile unit with connection for hot-air drying	•
Mobile unit for pipettes	•
Product characteristics	
Material	Stainless steel, Plastic
Colour	Stainless Steel, Black
Application area	
Reprocessing of pipettes / phials	•
Capacity	
Pipettes [number]	38
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in in. (mm)	19 13/16 (502)
External dimensions, net width in in. (mm)	21 1/8 (535)
External dimensions, net depth in in. (mm)	20 3/8 (517)
External dimensions, gross height in in. (mm)	19 11/16 (500)
External dimensions, gross width in in. (mm)	21 5/16 (540)
External dimensions, gross depth in in (mm)	20 7/8 (530)
Net weight in lb (kg)	14 (6.20)
Gross weight in lb (kg)	18 (8.10)