


E 103/1

Insert

For the optimum loading of test tubes up to 12 x 75 mm.



- Suitable for approx. 200 Test tubes (up to 12 x 75 mm)
- Divided into 6 surfaces
- The mesh width is 8 x 8 mm
- Including lid A 13

EAN: 4002514231051 / Material number: 06907630 / Old Material Number: 69510302D

Product category	
Glassware Washers, laboratory	•
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
G 7825	•
PG 8527	•
PG 8535	•
PG 8536	•
PG 8504	•
PG 8583	•
PG 8583 CD	•
PG 8593	•
PLW 8615	•
PLW 8616	•
PLW 8617	•
PLW 6011	•
PLW 6111	•
PLW 8505	•
Product category	
Insert / module for test tubes	•
Product characteristics	
Material	Stainless steel
Colour	Stainless Steel
Application area	
Reprocessing of test tubes	•
Reprocessing of laboratory glassware	•
Capacity	
Test tubes [number]	200
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in in. (mm)	4 13/16 (122)
External dimensions, net width in in. (mm)	7 7/8 (200)
External dimensions, net depth in in. (mm)	8 11/16 (220)
External dimensions, gross height in in. (mm)	5 9/16 (140)
External dimensions, gross width in in. (mm)	9 11/16 (245)
External dimensions, gross depth in in (mm)	9 7/8 (250)
Net weight in lb (kg)	2 (0.54)
Gross weight in lb (kg)	2 (0.74)
Accessory pack	
A 13 lids [number]	1