


E 104/1

Insert

For the optimum loading of test tubes up to 12 x 105 mm.



- Suitable for approx. 200 Test tubes (up to 12 x 105 mm)
- Capacity including lid A 13
- The mesh width is 8 x 8 mm

EAN: 4002514231068 / Material number: 06907640 / Old Material Number: 69510402D

<b>Product category</b>	
Glassware Washers, laboratory	•
Large-capacity washer-disinfectors, laboratory	•
<b>Machine affiliation</b>	
G 7825	•
PG 8527	•
PG 8535	•
PG 8536	•
PG 8504	•
PG 8583	•
PG 8583 CD	•
PG 8593	•
PLW 8615	•
PLW 8616	•
PLW 8617	•
PLW 6011	•
PLW 6111	•
PLW 8505	•
<b>Product category</b>	
Insert / module for test tubes	•
<b>Product characteristics</b>	
Material	Stainless steel
Colour	Stainless Steel
<b>Application area</b>	
Reprocessing of test tubes	•
Reprocessing of laboratory glassware	•
<b>Capacity</b>	
Test tubes [number]	200
<b>Standard electrical connection</b>	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
<b>Dimensions and weight</b>	
External dimensions, net height in in. (mm)	6 (152)
External dimensions, net width in in. (mm)	7 7/8 (200)
External dimensions, net depth in in. (mm)	9 1/16 (230)
External dimensions, gross height in in. (mm)	6 1/2 (165)
External dimensions, gross width in in. (mm)	9 1/2 (240)
External dimensions, gross depth in in (mm)	9 7/8 (250)
Net weight in lb (kg)	2 (0.60)
Gross weight in lb (kg)	2 (0.82)
<b>Accessory pack</b>	
A 13 lids [number]	1