


E 903/1

Module

for the optimum loading of TUR sets (transurethral resection).



- Holds up to 10 hollow instruments

EAN: 4002514256290 / Material number: 06943230 / Old Material Number: 69590302D

<b>Product category</b>	
Washer-disinfectors, medical	•
Large-capacity washer-disinfectors, medical	•
<b>Machine affiliation</b>	
G 7823	•
G 7824	•
G 7825	•
G 7826	•
PG 8527	•
PG 8528	•
PG 8536	•
<b>Product category</b>	
Insert / module for MIS instruments	•
<b>Product characteristics</b>	
Material	Stainless steel
Colour	Stainless Steel
<b>Application area</b>	
Reprocessing of MIS instruments	•
<b>Capacity</b>	
Lumened instruments [number]	10
<b>Standard electrical connection</b>	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
<b>Dimensions and weight</b>	
External dimensions, net height in in. (mm)	1 5/8 (40)
External dimensions, net width in in. (mm)	18 3/16 (461)
External dimensions, net depth in in. (mm)	20 1/8 (510)
External dimensions, gross height in in. (mm)	7 1/8 (180)
External dimensions, gross width in in. (mm)	21 11/16 (550)
External dimensions, gross depth in in (mm)	22 1/16 (560)
Net weight in lb (kg)	7 (2.93)
Gross weight in lb (kg)	9 (3.65)
<b>Accessory pack</b>	
E 442 irrigation sleeve, Ø 11 mm, length 121 mm, with cap Ø 6 mm [number]	3
E 444 drum for fibre optic cables / suction hoses [number]	1
E 447 female adapter, for male Luer Lock [number]	1
E 448 silicone hose with female Luer Lock adapter [number]	4
E 453 injector nozzle 4 x 110 mm with retaining clip [number]	3
E 454 injector nozzle for trocar sleeves 10–15 mm [number]	1
E 467 irrigation sleeve 205 mm for MIS instruments / clip forceps [number]	3
E 469 irrigation sleeve 300 mm for MIS instruments / urology [number]	3
E 907/1 mesh tray with lid for small parts [number]	1
2 m silicone hose, Ø 5 mm	•

E 903/1

Module

for the optimum loading of TUR sets (transurethral resection).

EAN: 4002514256290 / Material number: 06943230 / Old Material Number: 69590302D

Plastic support [number]

1

