

E 905/1

Module

for optimum loading of MIS instruments, such as arthroscopy and urology.



- Holds up to 16 hollow instruments
- Subdivision for arthroscopy and laparoscopy, for example.

EAN: 4002514256306 / Material number: 06943240 / Old Material Number: 69590502D

Product category	
Washer-disinfectors, medical	•
Large-capacity washer-disinfectors, medical	•
Machine affiliation	
G 7823	•
G 7824	•
G 7825	•
G 7826	•
PG 8527	•
PG 8528	•
PG 8536	•
Product category	
Insert / module for MIS instruments	•
Product characteristics	
Material	Stainless steel
Colour	Stainless Steel
Application area	
Reprocessing of MIS instruments	•
Capacity	
Lumened instruments [number]	16
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in in. (mm)	1 5/8 (40)
External dimensions, net width in in. (mm)	18 3/16 (461)
External dimensions, net depth in in. (mm)	20 1/8 (510)
External dimensions, gross height in in. (mm)	6 3/4 (170)
External dimensions, gross width in in. (mm)	21 5/16 (540)
External dimensions, gross depth in in (mm)	22 1/16 (560)
Net weight in lb (kg)	7 (3.01)
Gross weight in lb (kg)	8 (3.62)
Accessory pack	
E 336 MIBO irrigation sleeve, Ø 11 mm, length 121 mm [number]	1
E 362 blind stopper [number]	2
E 442 irrigation sleeve, Ø 11 mm, length 121 mm, with cap Ø 6 mm [number]	1
E 447 female adapter, for male Luer Lock [number]	2
E 448 silicone hose with female Luer Lock adapter [number]	4
E 449 male adapter, for female Luer Lock [number]	2
E 452 injector nozzle 2.5 x 60 mm [number]	4
E 453 injector nozzle 4 x 110 mm with retaining clip [number]	3
E 454 injector nozzle for trocar sleeves 10–15 mm [number]	3
E 464 holder for injector nozzle E 454 [number]	1
E 472 spring clip for injector nozzles Ø 4 mm [number]	1

E 905/1

Module

for optimum loading of MIS instruments, such as arthroscopy and urology.

EAN: 4002514256306 / Material number: 06943240 / Old Material Number: 69590502D



E 907/1 mesh tray with lid for small parts [number]

1