

E 906/1

Module

for optimum loading of MIS instruments, such as arthroscopy and urology.



- Holds up to 10 hollow instruments
- Subdivision for arthroscopy and laparoscopy, for example.

EAN: 4002514256313 / Material number: 06943250 / Old Material Number: 69590602D

<b>Product category</b>	
Washer-disinfectors, medical	•
Large-capacity washer-disinfectors, medical	•
<b>Machine affiliation</b>	
G 7823	•
G 7824	•
G 7825	•
G 7826	•
PG 8527	•
PG 8528	•
PG 8536	•
<b>Product category</b>	
Insert / module for MIS instruments	•
<b>Product characteristics</b>	
Material	Stainless steel
Colour	Stainless Steel
<b>Application area</b>	
Reprocessing of MIS instruments	•
<b>Capacity</b>	
Lumened instruments [number]	10
<b>Standard electrical connection</b>	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
<b>Dimensions and weight</b>	
External dimensions, net height in in. (mm)	1 5/8 (40)
External dimensions, net width in in. (mm)	18 3/16 (461)
External dimensions, net depth in in. (mm)	20 1/8 (510)
External dimensions, gross height in in. (mm)	6 5/16 (160)
External dimensions, gross width in in. (mm)	21 1/2 (545)
External dimensions, gross depth in in. (mm)	20 7/8 (530)
Net weight in lb (kg)	7 (2.87)
Gross weight in lb (kg)	8 (3.56)
<b>Accessory pack</b>	
E 336 MIBO irrigation sleeve, Ø 11 mm, length 121 mm [number]	1
E 362 blind stopper [number]	2
E 442 irrigation sleeve, Ø 11 mm, length 121 mm, with cap Ø 6 mm [number]	5
E 443 irrigation sleeve, Ø 11 mm, length 121 mm, with cap Ø 10 mm [number]	3
E 448 silicone hose with female Luer Lock adapter [number]	2
E 454 injector nozzle for trocar sleeves 10–15 mm [number]	1
E 456 spring stay for MIS instruments [number]	2
E 464 holder for injector nozzle E 454 [number]	1
E 908/1 insert for separable MIS instruments / work inserts [number]	1