

E 907/1

Mesh trays

for the optimum loading of different small-scale instruments.



- With lid for small parts
- Mesh width 3 x 1 mm
- Hook for hanging in E 905/1 and E 903/1

EAN: 4002514256320 / Material number: 06943820 / Old Material Number: 69590702D

Product category	
Washer-disinfectors, medical	•
Large-capacity washer-disinfectors, medical	•
Machine affiliation	
G 7823	•
G 7824	•
G 7825	•
G 7826	•
PG 8527	•
PG 8528	•
PG 8535	•
PG 8536	•
PG 8562	•
Product category	
Mesh tray	•
Product characteristics	
Material	Stainless steel
Colour	Stainless Steel
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in in. (mm)	1 15/16 (48)
External dimensions, net width in in. (mm)	5 1/8 (129)
External dimensions, net depth in in. (mm)	6 3/4 (170)
External dimensions, gross height in in. (mm)	2 13/16 (70)
External dimensions, gross width in in. (mm)	6 1/8 (155)
External dimensions, gross depth in in (mm)	9 1/16 (230)
Net weight in lb (kg)	1 (0.30)
Gross weight in lb (kg)	1 (0.37)