

E 912/2

Mobile unit

For optimum loading of 4 ISO containers including lids.



- Container size: 180 x 400 x 600 mm
- Magnetic strip for wash cart recognition (incl. magnets ML)

EAN: 4002514257976 / Material number: 06943860 / Old Material Number: 69591203D

Product category	
Large-capacity washer-disinfectors, medical	•
Machine affiliation	
PG 8527	•
PG 8528	•
Product category	
Mobile unit for containers	•
Product characteristics	
Material	Stainless steel
Colour	Stainless Steel
Application area	
Reprocessing of containers	•
Capacity	
ISO containers height 18 cm [quantity]	4
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in in. (mm)	17 15/16 (455)
External dimensions, net width in in. (mm)	25 1/4 (640)
External dimensions, net depth in in. (mm)	31 1/8 (790)
External dimensions, gross height in in. (mm)	25 5/8 (650)
External dimensions, gross width in in. (mm)	25 5/8 (650)
External dimensions, gross depth in in (mm)	31 1/2 (800)
Net weight in lb (kg)	29 (13.00)
Gross weight in lb (kg)	49 (22.00)
Accessory pack	
Magnetic strip for automatic program recognition (incl. magnets) [number]	1