

E 915/3

Mobile unit

for optimum loading of approx. 7 AN-sets or 20 breathing tubes and accessories.



- Holds breathing tubes on a spiral rack
- Injector system
- Equipped with connection for hot air drying
- Magnetic strip for mobile unit recognition (incl. magnets)

EAN: 4002514257983 / Material number: 06943870 / Old Material Number: 69591504D

Product category	
Large-capacity washer-disinfectors, medical	•
Machine affiliation	
PG 8527	•
PG 8528	•
Product category	
Mobile unit with connection for hot-air drying	•
Mobile unit for anaesthesia instruments	•
Product characteristics	
Material	Stainless steel
Colour	Stainless Steel
Application area	
Reprocessing anaesthetic instruments	•
Capacity	
AN sets [number]	7
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in in. (mm)	25 1/8 (637)
External dimensions, net width in in. (mm)	25 1/4 (640)
External dimensions, net depth in in. (mm)	31 1/8 (790)
External dimensions, gross height in in. (mm)	30 3/4 (780)
External dimensions, gross width in in. (mm)	26 7/16 (670)
External dimensions, gross depth in in (mm)	31 15/16 (810)
Net weight in lb (kg)	56 (25.00)
Gross weight in lb (kg)	56 (25.00)
Accessory pack	
E 142 insert 1/2, DIN mesh tray [number]	1
E 431 injector nozzle for concertina hose, Ø 8 x 193 mm [number]	1
E 466 injector nozzle 8 x 333 mm, for breathing bags [number]	10
E 496 injector nozzle 4 x 120 mm [number]	18
E 916 holder for 5 concertina breathing tubes [number]	1
E 917 holder for 5 silicone breathing tubes [number]	4
E 918 holder for 5 paediatric breathing tubes [number]	1
Spring holders, Ø 12 mm, for breathing tubes up to 1.5 m long [number]	20
Magnetic strip for automatic program recognition (incl. magnets) [number]	1