


## U 536 Basket

with low top frame for the optimum loading of 36 glasses up to 20 cm tall.




- With 36 compartments
- Compartment size 2 7/8 (73) in (mm) x 2 7/8 (73) in (mm) mm
- Includes extender set comprising short and long extenders
- Height adjustable top frame extendable up to 20 cm
- Width 50, depth 50 cm, width with Gastronorm (GN) sides 53 cm

EAN: 4002514278698 / Material number: 06984090 / Old Material Number: 67653601D

Product category	
Fresh-water dishwashers	•
Tank dishwashers	•
Throughfeed tank dishwashers	•
Machine affiliation	
G 8066	•
PG 8067	•
PG 8041	•
PG 8043	•
PG 8045	•
PG 8055	•
PG 8055 U	•
PG 8056	•
PG 8056 U	•
PG 8057 TD	•
PG 8057 TD U	•
PG 8058	•
PG 8058 U	•
PG 8059	•
PG 8059 U	•
PG 8060	•
PG 8061	•
PG 8061 U	•
PG 8165	•
PG 8166	•
PG 8169	•
PG 8172	•
PG 8172 Eco	•
PTD 901	•
PG 8096	•
PG 8096 U	•
PG 8099	•
PG 8099 U	•
PFD 400	•
PFD 400 U	•
PFD 401	•
PFD 401 U	•
PFD 402	•
PFD 402 U	•
PFD 404	•
PFD 404 U	•
PFD 405	•
PFD 405 U	•
PFD 407	•
PFD 407 U	•
PFD 408 U	•
PTD 702	•
PTD 703	•
PTD 704	•

U 536  
Basket

with low top frame for the optimum loading of 36 glasses up to 20 cm tall.



EAN: 4002514278698 / Material number: 06984090 / Old Material Number: 67653601D

Product category	
Basket for glasses	•
Product characteristics	
Material	Plastic
Colour	Blue, Green
Number of compartments	36
Compartment width in mm	2 7/8 (73)
Compartment length in mm	2 7/8 (73)
Capacity	
Glasses 110–140 mm high [number]	36
Glasses 175–200 mm high [number]	36
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in in. (mm)	9 3/16 (233)
External dimensions, net width in in. (mm)	19 11/16 (500)
External dimensions, net depth in in. (mm)	19 11/16 (500)
External dimensions, gross height in in. (mm)	10 1/4 (260)
External dimensions, gross width in in. (mm)	20 1/2 (520)
External dimensions, gross depth in in. (mm)	20 7/8 (530)
Net weight in lb (kg)	6 (2.70)
Gross weight in lb (kg)	8 (3.20)