

E 949

Module

for the optimum loading of 8 intramedullary nails and/or medullary reamers.



- Equipped with 8 holders for intramedullary nails*
- Equipped with 6 extension sleeves for medullary reamers**

* The holders are pre-installed.

** The extension sleeves are supplied loose.

EAN: 4002514312040 / Material number: 07025320 / Old Material Number: 69594901D

Product category	
Large-capacity washer-disinfectors, medical	•
Machine affiliation	
PG 8527	•
PG 8528	•
Product characteristics	
Material	Stainless steel
Colour	Stainless Steel
Application area	
Reprocessing of orthopaedic surgery instruments	•
Capacity	
Bone marrow drill [number]	8
Standard electrical connection	
Number of phases	0
Voltage in V	0-0
Frequency in Hz	0-0
Total rated load in kW	0.00-0.00
Fuse rating in A	0-0
Dimensions and weight	
External dimensions, net height in in. (mm)	4 9/16 (115)
External dimensions, net width in in. (mm)	13 7/8 (352)
External dimensions, net depth in in. (mm)	22 (558)
External dimensions, gross height in in. (mm)	11 7/16 (290)
External dimensions, gross width in in. (mm)	20 7/8 (530)
External dimensions, gross depth in in. (mm)	22 1/16 (560)
Net weight in lb (kg)	78 (35.00)
Gross weight in lb (kg)	100 (45.00)