


E 705/2

Mobile unit

for the optimum loading of MIS instruments up to 600 mm length.



- Up to 42 hollow instruments can be held
- Equipped with connection for hot air drying
- Magnetic strip for mobile unit recognition (without magnets ML)
- The vertical clearance for the 3 levels is different

EAN: 4002514313399 / Material number: 07029810 / Old Material Number: 69570503D

Product category	
Large-capacity washer-disinfectors, medical	•
Machine affiliation	
G 7826	•
Product category	
Mobile unit with connection for hot-air drying	•
Mobile unit for MIS instruments	•
Product characteristics	
Material	Stainless steel
Colour	Stainless Steel
Application area	
Reprocessing of MIS instruments	•
Capacity	
Lumened instruments [number]	42
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in in. (mm)	20 7/8 (530)
External dimensions, net width in in. (mm)	26 7/16 (670)
External dimensions, net depth in in. (mm)	23 5/8 (600)
External dimensions, gross height in in. (mm)	21 11/16 (550)
External dimensions, gross width in in. (mm)	25 1/4 (640)
External dimensions, gross depth in in. (mm)	25 5/8 (650)
Net weight in lb (kg)	28 (12.50)
Gross weight in lb (kg)	33 (14.60)
Accessory pack	
E 336 MIBO irrigation sleeve, Ø 11 mm, length 121 mm [number]	3
E 362 blind stopper [number]	2
E 442 irrigation sleeve, Ø 11 mm, length 121 mm, with cap Ø 6 mm [number]	15
E 443 irrigation sleeve, Ø 11 mm, length 121 mm, with cap Ø 10 mm [number]	5
E 444 drum for fibre optic cables / suction hoses [number]	1
E 445 – 12 irrigation sleeve covers E 442, opening with Ø 6 mm [number]	1
E 446 – 12 irrigation sleeve covers E 443, opening with Ø 10 mm [number]	1
E 447 female adapter, for male Luer Lock [number]	3
E 448 silicone hose with female Luer Lock adapter [number]	6
E 449 male adapter, for female Luer Lock [number]	5
E 451 insert 1/6 mesh tray [number]	1
E 452 injector nozzle 2.5 x 60 mm [number]	3
E 453 injector nozzle 4 x 110 mm with retaining clip [number]	8
E 454 injector nozzle for trocar sleeves 10–15 mm [number]	6

E 705/2

Mobile unit

for the optimum loading of MIS instruments up to 600 mm length.



EAN: 4002514313399 / Material number: 07029810 / Old Material Number: 69570503D

E 456 spring stay for MIS instruments [number]	4
E 464 holder for injector nozzle E 454 [number]	3
E 472 spring clip for injector nozzles Ø 4 mm [number]	2
2 m silicone hose, Ø 5 mm	•
Magnetic strip for automatic program recognition (without magnets) [number]	1