

TK 180
Tumble dryer basket
For hanging in the dryer.



- Suitable for PT 5186 and PT 7186
- For drying delicate fabrics
- Ideal for drying sports shoes for example

EAN: 4002514715698 / Material number: 07289550 / Old Material Number: 59018001D

Product category	
Tumble dryers	•
Machine / accessory category	
PT 5186	•
PT 7186	•
PT 7188	•
PWT 6089	•
Product category	
Dryer basket	•
Product characteristics	
Electrical connection	OHNE
Dimensions and weight	
External dimensions, net height in mm	0 (0)
External dimensions, net width in mm	0 (0)
External dimensions, net depth in mm	0 (0)
External dimensions, gross height in mm	18 3/4 (475)
External dimensions, gross width in mm	21 3/4 (550)
External dimensions, gross depth in mm	21 3/4 (550)
Net weight [kg]	3 (0.96)
Gross weight in kg	7 (2.94)