

E 478/1

Holder

for the optimum cleaning of 4 narrow-lumined cannulae (Sautter cannulae).



- Equipped with a stainless steel filter plate (FP-ED)

EAN: 4002514496948 / Material number: 07292920 / Old Material Number: 69747802D

Product classification	
Washer-disinfectors, medical	•
Large-capacity washer-disinfectors, medical	•
Machine affiliation	
G 7826	•
PG 8527	•
PG 8528	•
PG 8536	•
PG 8582	•
PG 8582 CD	•
PG 8592	•
Product category	
Holder	•
Product characteristics	
Material	Stainless steel
Colour	Stainless Steel
Application area	
Reprocessing of MIS instruments	•
Reprocessing of ophthalmology instruments	•
Capacity	
Lumened instruments [number]	4
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in mm	0
External dimensions, net width in mm	0
External dimensions, net depth in mm	0
External dimensions, gross height in mm	3 3/16 (80)
External dimensions, gross width in mm	3 3/4 (95)
External dimensions, gross depth in in (mm)	4 3/8 (110)
Net weight in kg	1 (0.19)
Gross weight in kg	1 (0.23)