

E 926
Module
for the optimum loading of long MIS instruments.



- Holds up to 13 hollow instruments

EAN: 4002514504353 / Material number: 07297050 / Old Material Number: 69592601D

Product category	
Large-capacity washer-disinfectors, medical	•
Machine affiliation	
PG 8527	•
PG 8528	•
Product category	
Insert / module for MIS instruments	•
Product characteristics	
Material	Stainless steel
Colour	Stainless Steel
Application area	
Reprocessing of MIS instruments	•
Capacity	
Lumened instruments [number]	13
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in in. (mm)	1 5/8 (40)
External dimensions, net width in in. (mm)	22 3/8 (568)
External dimensions, net depth in in. (mm)	20 1/8 (510)
External dimensions, gross height in in. (mm)	11 1/16 (280)
External dimensions, gross width in in. (mm)	22 7/8 (580)
External dimensions, gross depth in in (mm)	24 13/16 (630)
Net weight in lb (kg)	6 (2.40)
Gross weight in lb (kg)	6 (2.60)
Accessory pack	
E 336 MIBO irrigation sleeve, Ø 11 mm, length 121 mm [number]	3
E 362 blind stopper [number]	2
E 442 irrigation sleeve, Ø 11 mm, length 121 mm, with cap Ø 6 mm [number]	5
E 443 irrigation sleeve, Ø 11 mm, length 121 mm, with cap Ø 10 mm [number]	3
E 448 silicone hose with female Luer Lock adapter [number]	2
E 454 injector nozzle for trocar sleeves 10–15 mm [number]	1
E 456 spring stay for MIS instruments [number]	2
E 464 holder for injector nozzle E 454 [number]	1