


O 167/1

Upper basket

for the optimum loading of insoles for theatre shoes.



- Suitable for up to 28 insoles
- With integrated magnetic spray arm
- Magnet for use with spray arm sensors for PG 8535, PG 8536

EAN: 4002514556215 / Material number: 07333080 / Old Material Number: 69116703D

Product category	
Washer-disinfectors, medical	•
Machine affiliation	
PG 8535	•
PG 8536	•
Product category	
Insert for OR shoes	•
Upper basket	•
Product characteristics	
Material	Stainless steel
Colour	Stainless Steel
Application area	
Reprocessing of OR shoes	•
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in in. (mm)	7 11/16 (195)
External dimensions, net width in in. (mm)	18 15/16 (480)
External dimensions, net depth in in. (mm)	18 3/4 (475)
External dimensions, gross height in in. (mm)	8 5/16 (210)
External dimensions, gross width in in. (mm)	19 1/2 (495)
External dimensions, gross depth in in. (mm)	21 5/16 (540)
Net weight in lb (kg)	6 (2.40)
Gross weight in lb (kg)	8 (3.23)