

E 429/1

Mobile unit

for optimum cleaning of micro instruments (eye OP), for 2-4 OP sets.



- With integrated magnetic spray arm
- Level 1: e.g. 4 x insert E441/1 or 2 x insert E142 per level
- Level 2 left: 20 connections for hollow instruments
- Level 2 right: 16 Luer lock adapters, male, horizontal Hose
- Bracket for ML/2 for wash cart recognition / spray arm control

EAN: 4002514556253 / Material number: 07333120 / Old Material Number: 69542902D

Product category	
Washer-disinfectors, medical	•
Machine affiliation	
PG 8535	•
Product category	
Mobile unit for ophthalmology instruments	•
Product characteristics	
Material	Stainless steel
Colour	Stainless Steel
Application area	
Reprocessing of ophthalmology instruments	•
Capacity	
Surgery sets for micro instruments [number]	2-4
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in in. (mm)	20 1/2 (520)
External dimensions, net width in in. (mm)	21 5/16 (540)
External dimensions, net depth in in. (mm)	21 5/16 (540)
External dimensions, gross height in in. (mm)	20 7/8 (530)
External dimensions, gross width in in. (mm)	21 11/16 (550)
External dimensions, gross depth in in. (mm)	21 11/16 (550)
Net weight in lb (kg)	6 (2.40)
Gross weight in lb (kg)	8 (3.30)
Accessory pack	
E 476 holder [number]	1
E 477 stopper [number]	1
E 790 Luer Lock connection element female to female [number]	8
E 791 Luer Lock connection element female to male [number]	8