

E 440/3

Mobile unit

for optimum cleaning of micro instruments (eye OP), for 4 OP sets.



- With integrated magnetic spray arm
- Levels 1 + 2 e.g. 4 x E 441/1 or 2 x E 142 per level
- Level 3 left: 20 connections for hollow instruments
- Level 3 right: 16 Luer lock adapters, male, horizontal Hose
- Bracket for ML/2 for wash cart recognition / spray arm control

EAN: 4002514556260 / Material number: 07333140 / Old Material Number: 69544004D

Product category	
Washer-disinfectors, medical	•
Machine affiliation	
PG 8536	•
Product category	
Mobile unit for ophthalmology instruments	•
Product characteristics	
Material	Stainless steel
Colour	Stainless Steel
Application area	
Reprocessing of ophthalmology instruments	•
Capacity	
Surgery sets for micro instruments [number]	4
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in in. (mm)	19 13/16 (502)
External dimensions, net width in in. (mm)	21 1/8 (535)
External dimensions, net depth in in. (mm)	20 5/16 (515)
External dimensions, gross height in in. (mm)	21 5/16 (540)
External dimensions, gross width in in. (mm)	10 11/16 (270)
External dimensions, gross depth in in. (mm)	20 11/16 (525)
Net weight in lb (kg)	6 (2.40)
Gross weight in lb (kg)	8 (3.30)
Accessory pack	
E 476 holder [number]	1
E 477 stopper [number]	1
E 790 Luer Lock connection element female to female [number]	8
E 791 Luer Lock connection element female to male [number]	8