

E 942/3

Module

for the optimum loading of up to 116 pipettes.



- Suitable for 116 pipettes up to 580 mm in height
- Compartment size ?? x ??

EAN: 4002514605036 / Material number: 07459390 / Old Material Number: 69594204D

Product category	
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PG 8527	•
PG 8528	•
Product category	
Insert / module for pipettes	•
Product characteristics	
Material	Stainless steel
Colour	Stainless Steel
Application area	
Reprocessing of pipettes / phials	•
Capacity	
Pipettes [number]	116
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in in. (mm)	11 (279)
External dimensions, net width in in. (mm)	22 (558)
External dimensions, net depth in in. (mm)	13 7/8 (352)
External dimensions, gross height in in. (mm)	16 9/16 (420)
External dimensions, gross width in in. (mm)	21 11/16 (550)
External dimensions, gross depth in in (mm)	22 1/2 (570)
Net weight in lb (kg)	78 (35.00)
Gross weight in lb (kg)	100 (45.00)