

E 944/2

Module

for the optimum loading of narrow necked glassware 500 - 1,000 ml.



- 15 nozzles (E 352) 6 x 220 mm with clips (E 354)

EAN: 4002514605050 / Material number: 07459410 / Old Material Number: 69594403D

<b>Product category</b>	
Large-capacity washer-disinfectors, laboratory	•
<b>Machine affiliation</b>	
PG 8527	•
PG 8528	•
<b>Product category</b>	
Insert / module for narrow necked glassware	•
<b>Product characteristics</b>	
Material	Stainless steel
Colour	Stainless Steel
<b>Application area</b>	
Reprocessing of laboratory glassware	•
<b>Capacity</b>	
Narrow necked glassware [number]	15
<b>Standard electrical connection</b>	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
<b>Dimensions and weight</b>	
External dimensions, net height in in. (mm)	9 7/8 (250)
External dimensions, net width in in. (mm)	22 (558)
External dimensions, net depth in in. (mm)	13 7/8 (352)
External dimensions, gross height in in. (mm)	7 1/8 (180)
External dimensions, gross width in in. (mm)	21 11/16 (550)
External dimensions, gross depth in in (mm)	22 1/2 (570)
Net weight in lb (kg)	9 (3.92)
Gross weight in lb (kg)	10 (4.42)