


E 947/2

Module

for cleaning 88 centrifuge tubes, vials or test tubes, for example.



- Equipped with 88 nozzles, 2.5 x 110 mm

EAN: 4002514605074 / Material number: 07459430 / Old Material Number: 69594703D

Product category	
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PG 8527	•
PG 8528	•
Product category	
Insert / module for test tubes	•
Product characteristics	
Material	Stainless steel, Plastic
Colour	Stainless Steel, Black
Application area	
Reprocessing of pipettes / phials	•
Capacity	
Test tubes [number]	88
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in in. (mm)	6 3/4 (170)
External dimensions, net width in in. (mm)	22 (558)
External dimensions, net depth in in. (mm)	13 7/8 (352)
External dimensions, gross height in in. (mm)	11 7/16 (290)
External dimensions, gross width in in. (mm)	20 7/8 (530)
External dimensions, gross depth in in. (mm)	22 1/16 (560)
Net weight in lb (kg)	12 (5.00)
Gross weight in lb (kg)	14 (6.00)