

E 555/2

Mobile unit

for cleaning instruments on 2 levels.



- Loading volume from below: H 210, W 487, D 592 mm
- Loading volume from above: H 219, W 495, D 548 mm
- Magnetic strip for mobile unit recognition (without magnets ML)
- Wide range of fitments available

EAN: 4002514606248 / Material number: 07460290 / Old Material Number: 69555503D

<b>Product category</b>	
Large-capacity washer-disinfectors, medical	•
<b>Machine affiliation</b>	
G 7823	•
G 7824	•
<b>Product category</b>	
Mobile unit for inserts	•
<b>Product characteristics</b>	
Material	Stainless steel
Colour	Stainless Steel
<b>Application area</b>	
Reprocessing of basins	•
Reprocessing of kidney dishes	•
Reprocessing of baby bottles	•
Reprocessing of surgical instruments	•
Reprocessing of containers	•
<b>Capacity</b>	
DIN mesh trays [number]	4
SPRI mesh trays [number]	2
<b>Standard electrical connection</b>	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
<b>Dimensions and weight</b>	
External dimensions, net height in in. (mm)	11 3/4 (298)
External dimensions, net width in in. (mm)	20 7/8 (530)
External dimensions, net depth in in. (mm)	25 7/16 (645)
External dimensions, gross height in in. (mm)	25 5/8 (650)
External dimensions, gross width in in. (mm)	21 11/16 (550)
External dimensions, gross depth in in (mm)	25 1/4 (640)
Net weight in lb (kg)	23 (10.00)
Gross weight in lb (kg)	27 (12.00)
<b>Accessory pack</b>	
Magnetic strip for automatic program recognition (without magnets) [number]	1