

E 521/2

Insert

for the optimum loading of 7 extraction, orthodontic and/or technical forceps.



- Compartment size: 21 x 80 mm
- Also suitable for forceps, such as those used in podiatry

EAN: 4002514616438 / Material number: 07476430 / Old Material Number: 68552103D

Product category	
Thermal disinfectors	•
Washer-disinfectors, medical	•
Machine affiliation	
G 7831	•
PG 8581	•
PG 8591	•
PWD 8531	•
PWD 8532	•
PWD 8534	•
Product category	
Insert / module for dental instruments	•
Product characteristics	
Material	Stainless steel
Colour	Stainless Steel
Application area	
Reprocessing of dental instruments	•
Capacity	
Forceps [number]	7
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in mm	5 3/8 (135) in (mm)
External dimensions, net width in mm	3 15/16 (100) in (mm)
External dimensions, net depth in mm	7 1/2 (189) in (mm)
External dimensions, gross height in mm	4 3/4 (120)
External dimensions, gross width in mm	6 15/16 (175)
External dimensions, gross depth in in (mm)	8 11/16 (220)
Net weight in kg	1 (0.42)
Gross weight in kg	2 (0.46)