

DOS NA 120

Bellows dispenser pump

For dispensing liquid neutral/alkaline agents.



- For dispensing liquid cleaning detergent, emulgator
- For dispensing chemical disinfectant
- For retro-fitting by Miele Service

EAN: 4002514641911 / Material number: 07514610 / Old Material Number: 69226005D

Product classification	
Glassware Washers, laboratory	•
Washer-disinfectors, medical	•
Machine affiliation	
PG 8536	•
Product category	
Dispenser system	•
Standard electrical connection	
Number of phases	0
Voltage in V	0
Frequency in Hz	0
Total rated load in kW	0.00
Fuse rating in A	0
Dimensions and weight	
External dimensions, net height in mm	0
External dimensions, net width in mm	0
External dimensions, net depth in mm	0
External dimensions, gross height in mm	11 13/16 (300)
External dimensions, gross width in mm	9 7/8 (250)
External dimensions, gross depth in in (mm)	16 9/16 (420)
Net weight in kg	6 (2.45)
Gross weight in kg	7 (2.89)